



Anti-ZBTB34 polyclonal antibody (DPAB-DC1147)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZBTB34 (zinc finger and BTB domain containing 34) is a protein-coding gene. GO annotations related to this gene include DNA binding and zinc ion binding. An important paralog of this gene is ZBTB22.
Specificity	Specific to recombinant protein GX1662. This antibody detects mZBTB34 protein. It also recognizes human ZBTB34 protein.
Immunogen	Recombinant GST fusion protein corresponding to 205 amino acids of mouse Zbtb34.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Zbtb34 zinc finger and BTB domain containing 34 [Mus musculus (house mouse)]
Official Symbol	ZBTB34
Synonyms	ZBTB34; zinc finger and BTB domain containing 34; Gm349; mKIAA1993; zinc finger and BTB domain-containing protein 34;
Entrez Gene ID	241311
Protein Refseq	NP_001078976
UniProt ID	A2ATY4
Chromosome Location	2 B; 2
Function	metal ion binding; molecular_function;