



Anti-CNKSR2 polyclonal antibody (DPAB-DC1155)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CNKSR2 (connector enhancer of kinase suppressor of Ras 2) is a protein-coding gene. Diseases associated with CNKSR2 include mastitis, and neuritis. GO annotations related to this gene include phospholipid binding. An important paralog of this gene is IPCEF1.
Specificity	Specific to recombinant protein GX0244. This antibody detects mCNKSR2 protein.
Immunogen	Recombinant GST fusion protein corresponding to 171 mouse Cnksr2. The sequence is VSACDPQDDIQPPEVEEEEEEEEEEAAGENVGEKNENREEKLGDSLQDLYRALEEASLSP LGHRISTKMEYKLSFIKRCNDPVMNEKLHRLRLKSTLKAREGEVAIIDKVLDNPD LTS KEFQQWKQMYLDLFLDICQSTTSNDPLSISSEVDVLTSSLTHTHSYIETHV
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Cnksr2 connector enhancer of kinase suppressor of Ras 2 [Mus musculus (house mouse)]
Official Symbol	CNKSR2
Synonyms	CNKSR2; connector enhancer of kinase suppressor of Ras 2; Cnk2; Ksr2; mKIAA0902; C030004A01; connector enhancer of kinase suppressor of ras 2; CNK homolog protein 2; connector enhancer of KSR2; connector enhancer of KSR 2;
Entrez Gene ID	245684
Protein Refseq	NP_808419
UniProt ID	Q80YA9
Chromosome Location	X F4; X
Function	identical protein binding; protein kinase binding;