



Anti-GAS2L2 (aa 665-776) polyclonal antibody (DPAB-DC1159)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene appears to crosslink microtubules and microfilaments and may be part of the cytoskeleton. This gene is mainly expressed in skeletal muscle.
Immunogen	GAS2L2 (NP_644814, 665 a.a. ~ 776 a.a) partial recombinant protein with GST tag. The sequence is LLKVDLEAWKAAPTGSPKPAVTPGPGSLKGKLGARQSGPRTKASLSAKGTHMRKVPPQGG QDCSASTVSASPEAPTPSPLDPNSDKAKACLSKGRRTLRLKPKRVPSIYKLLK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GAS2L2 growth arrest-specific 2 like 2 [Homo sapiens (human)]
Official Symbol	GAS2L2

Synonyms	GAS2L2; growth arrest-specific 2 like 2; GAR17; GAS2-like protein 2; GAS2-related protein on chromosome 17;
Entrez Gene ID	246176
Protein Refseq	NP_644814
UniProt ID	Q8NHY3
Chromosome Location	17q12
Function	actin filament binding; cytoskeletal adaptor activity; microtubule binding;