



Anti-ADM polyclonal antibody (DPAB-DC1169)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	involved in the regulation of adrenal development; plays a role in inhibition of apoptosis [RGD, Feb 2006]
Immunogen	A synthetic peptide (conjugated with KLH) corresponding to rat Adm. The sequence is YRQSMNQGSRSTGCRFGTCTMQKLAHQIYQFTDKDKDGMAPRNKISPQGY
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Adm adrenomedullin [Rattus norvegicus (Norway rat)]
Official Symbol	ADM

Synonyms	ADM; adrenomedullin; Ap; RATAP; H39316;
Entrez Gene ID	25026
Protein Refseq	NP_036847
UniProt ID	P43145
Chromosome Location	1q33
Pathway	Calcitonin-like ligand receptors; G alpha (s) signalling events; GPCR ligand binding; Signal Transduction
Function	adrenomedullin receptor binding; hormone activity; protein binding;