



Anti-PRP2 (aa 499-517) polyclonal antibody (DPAB-DC1190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 499-517 of <i>Saccharomyces cerevisiae</i> mis11. The sequence is RTIVIAFYGEDCYKANAW
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Yeast
Conjugate	Unconjugated
Applications	WB, IP,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	prp2 U2AF large subunit (U2AF-59) [Schizosaccharomyces pombe]
Official Symbol	PRP2
Synonyms	PRP2; mis11; U2AF large subunit (U2AF-59);

Entrez Gene ID	2539645
Protein Refseq	NP_595396
UniProt ID	P36629
Chromosome Location	chromosome: II
Pathway	Spliceosome.
Function	mRNA binding; snoRNA binding;