



Anti-ULP1 (aa 7-20) polyclonal antibody (DPAB-DC1194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 7-20 of Schizosaccharomyces pombe Ulp1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Schizosaccharomyces pombe
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ulp1 SUMO deconjugating enzyme Ulp1 [Schizosaccharomyces pombe]
Official Symbol	ULP1
Synonyms	ULP1; SUMO deconjugating enzyme Ulp1;

Entrez Gene ID	2540734
Protein Refseq	NP_595975
UniProt ID	O42957
Chromosome Location	chromosome: II
Function	SUMO-specific protease activity;