



Anti-SGO1 polyclonal antibody (DPAB-DC1195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to Schizosaccharomyces pombe Sgo1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pombe
Conjugate	Unconjugated
Applications	WB (Cell lysate), ICC, IF,
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	sgo1 inner centromere protein, shugoshin, Sgo1 [Schizosaccharomyces pombe]
Official Symbol	SGO1
Synonyms	SGO1; inner centromere protein, shugoshin, Sgo1;
Entrez Gene ID	2541326

Protein Refseq	NP_595378
UniProt ID	Q9P7A0
Chromosome Location	chromosome: II
Function	molecular_function;