



Anti-SWI6 (aa 314-328) polyclonal antibody (DPAB-DC1198)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 314-328 of Schizosaccharomyces pombe swi6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pombe
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	swi6 chromodomain protein Swi6 [Schizosaccharomyces pombe]
Official Symbol	SWI6

Synonyms	SWI6; SPAC824.10c; chromodomain protein Swi6;
Entrez Gene ID	2541633
Protein Refseq	NP_593449
UniProt ID	P40381
Chromosome Location	chromosome: I
Function	RNA binding; chromatin binding; methylated histone binding; nucleosomal histone binding