



Anti-RNF169 (aa 643-742) polyclonal antibody (DPAB-DC1201)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF169 (ring finger protein 169) is a protein-coding gene. GO annotations related to this gene include nucleosome binding and ligase activity.
Immunogen	KIAA1991 (XP_495886, 643 a.a. ~ 742 a.a) partial recombinant protein with GST tag. The sequence is SLSEEPLPSLRRGRKRHCKTKHLEQNGSLKKLRQTSGEVGLAPTD PVLREMEQKLQEEEE DRQLALQLQRMFDNERRTVSRRKGSVDQYLLRSSNMAGAK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNF169 ring finger protein 169 [Homo sapiens (human)]
Official Symbol	RNF169

Synonyms	RNF169; ring finger protein 169; E3 ubiquitin-protein ligase RNF169;
Entrez Gene ID	254225
Protein Refseq	NP_001092108
UniProt ID	Q8NCN4
Chromosome Location	11q13.4
Function	K63-linked polyubiquitin binding; ligase activity; nucleosome binding; zinc ion binding