



Anti-ULP2 (aa 622-635) polyclonal antibody (DPAB-DC1202)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 622-635 of Schizosaccharomyces pombe Ulp2. The sequence is SNNERQSLSSGSND
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Schizosaccharomyces pombe
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ulp2 SUMO deconjugating cysteine peptidase Ulp2 (predicted) [Schizosaccharomyces pombe]
Official Symbol	ULP2

Synonyms	ULP2; SUMO deconjugating cysteine peptidase Ulp2 (predicted);
Entrez Gene ID	2542268
Protein Refseq	NP_593475
UniProt ID	O13769
Chromosome Location	chromosome: I
Function	SUMO-specific protease activity;