



Anti-RAD22 (full length) polyclonal antibody (DPAB-DC1203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Schizosaccharomyces pombe</i> rad22.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Schizosaccharomyces pombe</i>
Conjugate	Unconjugated
Applications	WB, IF, IP,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.09% sodium azide
Preservative	0.09% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	rad22 DNA recombination protein Rad22 [<i>Schizosaccharomyces pombe</i>]
Official Symbol	RAD22
Synonyms	RAD22; rad22A; DNA recombination protein Rad22;

Entrez Gene ID	2543086
Protein Refseq	NP_593207
UniProt ID	P36592
Chromosome Location	chromosome: I
Pathway	Assembly of the RAD51-ssDNA nucleoprotein complex; Double-Strand Break Repair; Homologous Recombination Repair; Homologous recombination.
Function	damaged DNA binding; protein binding; single-stranded DNA binding;