



Anti-RAD4 (N-terminal) polyclonal antibody (DPAB-DC1205)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with <i>S. pombe</i> Cut5.
Immunogen	Recombinant GST fusion protein corresponding to N-terminus half <i>Schizosaccharomyces pombe</i> Cut5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. pombe</i>
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [rad4 BRCT domain protein Rad4 \[*Schizosaccharomyces pombe* \]](#)

Official Symbol	RAD4
Synonyms	RAD4; cut5; dpb11; dre3; BRCT domain protein Rad4;
Entrez Gene ID	2543281
Protein Refseq	NP_593190
UniProt ID	P32372
Chromosome Location	chromosome: I
Function	protein binding; protein complex scaffold;