



Anti-RHP51 (full length) polyclonal antibody (DPAB-DC1206)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length <i>Schizosaccharomyces pombe</i> rad51.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Schizosaccharomyces pombe</i>
Conjugate	Unconjugated
Applications	WB, IF, IP,
Format	Liquid
Size	50 µl
Buffer	In buffer containing 0.09% sodium azide
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	rhp51 RecA family recombinase Rhp51 [<i>Schizosaccharomyces pombe</i>]
Official Symbol	RHP51
Synonyms	RHP51; rad51; RecA family recombinase Rhp51;

Entrez Gene ID	2543580
Protein Refseq	NP_593882
Chromosome Location	chromosome: I
Pathway	Homologous recombination.
Function	ATP binding; ATPase activity; damaged DNA binding; double-stranded DNA binding