



Anti-LOC257054 (full length) polyclonal antibody (DPAB-DC1231)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MGC25181 (AAH31817, 1 a.a. ~ 220 a.a) full-length recombinant protein with GST tag. The sequence is MLGEILSAFEFMDAVCMQLVGRHLHLASPVQESPFYVLIETSGSNAGHDAEKLGHFLEHA LGSGLVTDGTMATDQRKVKMLWALRERITEALSRDGYVYKYDLSLPVERLYDIVTDLRAR LGPHAKHVVGYGHLGDGNLHLNVTAEAFSPSLLAALEPHVYEWTAGQQGSVSAEHGVGFR KRDVLGYSKPPGALQLMQQLKALLDPKGILNPYKTLPSQA
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	LOC257054 similar to D-2-hydroxyglutarate dehydrogenase [Homo sapiens (human)]
Official Symbol	LOC257054

Synonyms	LOC257054; similar to D-2-hydroxyglutarate dehydrogenase;
Entrez Gene ID	257054
Chromosome Location	2q37.3