



Anti-RNF214 (aa 541-640) polyclonal antibody (DPAB-DC1234)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF214 (ring finger protein 214) is a protein-coding gene. Diseases associated with RNF214 include migraine. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is TTC3.
Immunogen	DKFZp547C195 (NP_997226, 541 a.a. ~ 640 a.a) partial recombinant protein with GST tag. The sequence is LAEHERVAASTQPLGRIRALFPAPLAQISTPMFLPSAQVSYSGRSSHAPATCKLCLMCQK LVQPSELHPMACTHVLHKECIKFWAQTNTNDTCPCPTLK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RNF214 ring finger protein 214 \[Homo sapiens \(human\) \]](#)

Official Symbol	RNF214
Synonyms	RNF214; ring finger protein 214; RING finger protein 214;
Entrez Gene ID	257160
Protein Refseq	NP_001070707
UniProt ID	A0A024R3D4
Chromosome Location	11q23.3
Function	zinc ion binding;