



Anti-CCNE1 (C-terminal) polyclonal antibody (DPAB-DC1237)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	a cyclin-like protein; may be involved in growth stimulation [RGD, Feb 2006]
Immunogen	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of rat Ccne1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	200 µg
Buffer	In PBS (0.2% gelatin, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. Do not freeze.

GENE INFORMATION

Gene Name	Ccne1 cyclin E1 [Rattus norvegicus (Norway rat)]
Official Symbol	CCNE1

Synonyms	CCNE1; cyclin E1; Ccne; CYCLE; G1/S-specific cyclin-E1; cyclin E;
Entrez Gene ID	25729
Protein Refseq	NP_001094291
UniProt ID	B1WC54
Chromosome Location	1q22
Pathway	ATM; B Cell Receptor Signaling Pathway; Cell Cycle Checkpoints; Cell cycle
Function	cyclin-dependent protein serine/threonine kinase regulator activity; kinase activity; protein complex binding; protein kinase binding