



Anti-HECTD1 (aa 3-110) polyclonal antibody (DPAB-DC1245)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HECTD1 (HECT domain containing E3 ubiquitin protein ligase 1) is a protein-coding gene. GO annotations related to this gene include ubiquitin-protein ligase activity and metal ion binding. An important paralog of this gene is TRIP12.
Immunogen	HECTD1 (NP_056197, 3 a.a. ~ 110 a.a) partial recombinant protein with GST tag. The sequence is DVPDPTLLEWLQMGQGDERDMQLIALEQLCMLLLMSDNVDRCFETCPPRTFLPALCKIFL DESAPDNVLEVTARAITYYLDVSAECTRRIVGV DGAIKALCNRLVVVE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [HECTD1 HECT domain containing E3 ubiquitin protein ligase 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	HECTD1
Synonyms	HECTD1; HECT domain containing E3 ubiquitin protein ligase 1; EULIR; E3 ubiquitin-protein ligase HECTD1; E3 ligase for inhibin receptor; HECT domain-containing protein 1;
Entrez Gene ID	25831
Protein Refseq	NP_056197
UniProt ID	Q9ULT8
Chromosome Location	14q12
Function	ligase activity; metal ion binding; protein binding; ubiquitin-protein transferase activity