



Anti-ZNF345 (aa 1-111) polyclonal antibody (DPAB-DC1249)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNF345 (zinc finger protein 345) is a protein-coding gene. GO annotations related to this gene include DNA binding and zinc ion binding. An important paralog of this gene is ZNF79.
Immunogen	ZNF345 (NP_003410, 1 a.a. ~ 111 a.a) partial recombinant protein with GST tag. The sequence is MENLTKHSIECSSFRGDWECKNQFERKQGSQEGHFSEMIFTPEDMPTFSIQHQRIHTDEK LLECKECKGKDFSFVSVLVRHQRIHTGEKPYECKCGKAFGSGANLAYHQRI
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ZNF345 zinc finger protein 345 [Homo sapiens (human)]
Official Symbol	ZNF345

Synonyms	ZNF345; zinc finger protein 345; HZF10; zinc finger protein HZF10;
Entrez Gene ID	25850
Protein Refseq	NP_001229401
UniProt ID	B2RCE4
Chromosome Location	19q13.12
Function	DNA binding; metal ion binding; sequence-specific DNA binding transcription factor activity;