



Anti-RNF167 (aa 78-172) polyclonal antibody (DPAB-DC1267)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF167 is an E3 ubiquitin ligase that interacts with TSSC5 (SLC22A18; MIM 602631) and, together with UBCH6 (UBE2E1; MIM 602916), facilitates TSSC5 polyubiquitylation (Yamada and Gorbsky, 2006 [PubMed 16314844]).[supplied by OMIM, Mar 2008]
Immunogen	RNF167 (NP_056343, 78 a.a. ~ 172 a.a) partial recombinant protein with GST tag. The sequence is VNGSVFIALLRRFDCNFDLKVLAQKAGYGAAVVHNVNSNELLNMVWNSEEIQQQIWIPS VFIGERSSEYLRALFVYEKGARVLLVPDNTFPLGY
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RNF167 ring finger protein 167 \[Homo sapiens \(human\) \]](#)

Official Symbol	RNF167
Synonyms	RNF167; ring finger protein 167; LP2254; RING105; 5730408C10Rik; E3 ubiquitin-protein ligase RNF167;
Entrez Gene ID	26001
Protein Refseq	NP_056343
UniProt ID	Q9H6Y7
Chromosome Location	17p13.2
Function	ligase activity; protein binding; ubiquitin-protein transferase activity; zinc ion binding