



Anti-SUSD5 (aa 284-383) polyclonal antibody (DPAB-DC1269)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SUSD5 (sushi domain containing 5) is a protein-coding gene. GO annotations related to this gene include hyaluronic acid binding.
Immunogen	SUSD5 (XP_171054, 284 a.a. ~ 383 a.a) partial recombinant protein with GST tag. The sequence is TTVCCKSGSEQQIMRAVDVRIESNPVPGGTYSALCIKDEEKPCGDPPSPHTILQGRGTGL EMGDELLYVCAPGHIMGHRETAFTLLCNSCGEWYGLVQAC
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SUSD5 sushi domain containing 5 [Homo sapiens (human)]
Official Symbol	SUSD5

Synonyms	SUSD5; sushi domain containing 5; sushi domain-containing protein 5;
Entrez Gene ID	26032
Protein Refseq	NP_056366
UniProt ID	O60279
Chromosome Location	3p22.3
Function	hyaluronic acid binding;