



Anti-SHPRH (aa 515-535) polyclonal antibody (DPAB-DC1325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SHPRH (SNF2 histone linker PHD RING helicase, E3 ubiquitin protein ligase) is a protein-coding gene. GO annotations related to this gene include helicase activity and ligase activity.
Immunogen	A synthetic peptide corresponding to amino acids 515-535 of the mouse Shprh. The sequence is DKTKKQAVGSPRKIEKELRKS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [Shprh SNF2 histone linker PHD RING helicase \[Mus musculus \(house mouse\) \]](#)

Official Symbol	SHPRH
Synonyms	SHPRH; SNF2 histone linker PHD RING helicase; AA450458; AU024614; BC006883; E130018M05; 2610103K11Rik; D230017O13Rik; E3 ubiquitin-protein ligase SHPRH; SNF2, histone-linker, PHD and RING finger domain-containing helicase;
Entrez Gene ID	268281
Protein Refseq	NP_001071175
UniProt ID	Q7TPQ3
Chromosome Location	10 A1; 10
Function	ATP binding; DNA binding; UDP-N-acetylmuramoylalanyl-D-glutamyl-2,6-diaminopimelate-D-alanyl-D-alanine ligase activity; coenzyme F420-0 gamma-glutamyl ligase activity