



Anti-Glucose Oxidase (native) polyclonal antibody (DPAB-DC1326)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified <i>Aspergillus nige</i> Glucose Oxidase.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Fungi
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	MGG_07580 glucose oxidase [Magnaporthe oryzae 70-15]
Official Symbol	MGG_07580

Synonyms	MGG_07580; glucose oxidase;
Entrez Gene ID	2683500
UniProt ID	G4N2I1
Chromosome Location	chromosome: 3