



Anti-TETR (aa 31-49) polyclonal antibody (DPAB-DC1360)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 31-49 of Escherichia coli Tet Repressor. The sequence is AQKLGVEQPTLYWHVKNKR
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Escherichia coli
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	tetR hypothetical protein [Escherichia coli]
Official Symbol	TETR

Synonyms	TETR; hypothetical protein;
Entrez Gene ID	2716476
UniProt ID	P03038
Chromosome Location	plasmid: pC15-1a