



Anti-ADAMDEC1 (aa 361-470) polyclonal antibody (DPAB-DC1369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells.
Immunogen	ADAMDEC1 (NP_055294, 361 a.a. ~ 470 a.a) partial recombinant protein with GST tag. The sequence is PDVPFNTKCPSGSCVMNQYLSSKFPKDFSTSCRAHFERYLLSQKPKCLLQAPIPTNIMTT PVCGNHLLLEVGEDCDCGSPKECTNLCCEALTCKLKPGTDCGGDAPNHTTE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ADAMDEC1 ADAM-like, decysin 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	ADAMDEC1
Synonyms	ADAMDEC1; ADAM-like, decysin 1; M12.219; ADAM DEC1; decysin; disintegrin protease; ADAM-like protein decysin-1; a disintegrin and metalloproteinase domain-like protein decysin-1;
Entrez Gene ID	27299
Protein Refseq	NP_001138743
UniProt ID	B7Z6V5
Chromosome Location	8p21.2
Function	metalloendopeptidase activity; zinc ion binding;