



Anti-UTP20 (aa 946-1053) polyclonal antibody (DPAB-DC1377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	UTP20 is a component of the U3 small nucleolar RNA (snoRNA) (SNORD3A; MIM 180710) protein complex (U3 snoRNP) and is involved in 18S rRNA processing (Wang et al., 2007 [PubMed 17498821]).[supplied by OMIM, Jun 2009]
Immunogen	DRIM (NP_055318, 946 a.a. ~ 1053 a.a) partial recombinant protein with GST tag. The sequence is LHQDQMVQKITLDCIMTYKHPHVLPHYRENLRLLLED RSFKEEIVHFSISEDNAVVKTAHR ADLFPILMRILYGRMKNKTGSKTQGKSASGTRMAIVLRFLAGTQPEEI
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	UTP20 UTP20, small subunit (SSU) processome component, homolog (yeast) [Homo sapiens (human)]
------------------	---

Official Symbol	UTP20
Synonyms	UTP20; UTP20, small subunit (SSU) processome component, homolog (yeast); DRIM; small subunit processome component 20 homolog; NNP73; novel nucleolar protein 73; down regulated in metastasis; down-regulated in metastasis protein;
Entrez Gene ID	27340
Protein Refseq	NP_055318
UniProt ID	O75691
Chromosome Location	12q23
Function	poly(A) RNA binding; protein binding;