



Anti-MKLN1 (N-terminal) polyclonal antibody (DPAB-DC1383)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	MKLN1 (muskelin 1, intracellular mediator containing kelch motifs) is a protein-coding gene. Diseases associated with MKLN1 include ocular albinism, and albinism.
Immunogen	A synthetic peptide corresponding to N-terminus of mouse Mkl1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Mkl1 muskelin 1, intracellular mediator containing kelch motifs [Mus musculus (house mouse)]
------------------	--

Official Symbol	MKLN1
Synonyms	MKLN1; muskelin 1, intracellular mediator containing kelch motifs; AU015903; BB154892; A130067F06Rik; muskelin;
Entrez Gene ID	27418
Protein Refseq	NP_038819
UniProt ID	O89050
Chromosome Location	6; 6 B1-B2
Function	protein binding;