



Anti-GML (aa 48-158) polyclonal antibody (DPAB-DC1390)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GML (glycosylphosphatidylinositol anchored molecule like) is a protein-coding gene. Diseases associated with GML include ainhum, and dental pulp necrosis.
Immunogen	GML (NP_002057, 48 a.a. ~ 158 a.a) partial recombinant protein with GST tag. The sequence is CPYHIRRCMTISIRINSRELLVYKNCTNNCTFVYAAEQPPEAPGKIFKTNSFYWVCCNS MVCNAGGPTNLERDMLPDEVTEELPEGTVRLGVSKLLLSFASIIVSNILP
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GML glycosylphosphatidylinositol anchored molecule like [Homo sapiens (human)]
Official Symbol	GML

Synonyms	GML; glycosylphosphatidylinositol anchored molecule like; LY6DL; glycosyl-phosphatidylinositol-anchored molecule-like protein; GPI anchored molecule like protein; Glycosylphosphatidylinositol-anchored molecule-like protein; glycosylphosphatidylinositol anchored molecule like protein;
Entrez Gene ID	2765
Protein Refseq	NP_002057
UniProt ID	Q99445
Chromosome Location	8q24.3