



Anti-FAM3C (N-terminal) polyclonal antibody (DPAB-DC1397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM3C (family with sequence similarity 3, member C) is a protein-coding gene. Diseases associated with FAM3C include osteopoikilosis, and coffin-siris syndrome. GO annotations related to this gene include cytokine activity. An important paralog of this gene is FAM3D.
Immunogen	A synthetic peptide corresponding to N-terminus of mouse Fam3c.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Fam3c family with sequence similarity 3, member C [Mus musculus (house mouse)]
------------------	--

Official Symbol	FAM3C
Synonyms	FAM3C; family with sequence similarity 3, member C; Ilei; D6Wsu176e; protein FAM3C; interleukin-like EMT inducer;
Entrez Gene ID	27999
Protein Refseq	NP_613053
UniProt ID	Q91VU0
Chromosome Location	6 A3.1; 6 9.24 cM
Function	molecular_function;