



Anti-LOC283155 (aa 207-297) polyclonal antibody (DPAB-DC1418)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LOC283155 (XP_208545, 207 a.a. ~ 297 a.a) partial recombinant protein with GST tag. The sequence is RRCITGGIVEEYQLPYYNMVPSPSYEDMREVVCVKRLRPIVSNRWNSDECLRAVLKLMS ECWAHNPASRLTALRIKKMLAKMVESQDVKI
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	LOC283155 similar to ALK-3 [Homo sapiens (human)]
Official Symbol	LOC283155
Synonyms	LOC283155;
Entrez Gene ID	283155

