



Anti-HECTD4 (aa 3105-3209) polyclonal antibody (DPAB-DC1421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HECTD4 (HECT domain containing E3 ubiquitin protein ligase 4) is a protein-coding gene. Diseases associated with HECTD4 include esophageal squamous cell carcinoma, and esophagitis. GO annotations related to this gene include ubiquitin-protein ligase activity.
Immunogen	C12orf51 (XP_497354, 3105 a.a. ~ 3209 a.a) partial recombinant protein with GST tag. The sequence is DQHIEFFWGALEMFTQEELCKFIKFACNQERIPFTCPCKDGGPDTAHVPPYPMKIAPPDG TAGSPDSRYIRVETCMFMIKLPQYSSLEIMLEKLRCIAIHYREDPL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	HECTD4 HECT domain containing E3 ubiquitin protein ligase 4 [Homo sapiens (human)]
------------------	--

Official Symbol	HECTD4
Synonyms	HECTD4; HECT domain containing E3 ubiquitin protein ligase 4; POTAGE; C12orf51; probable E3 ubiquitin-protein ligase HECTD4; transmembrane protein C12orf51; HECT domain-containing protein 4; AF-1 specific protein phosphatase; probable E3 ubiquitin-protein ligase C12orf51;
Entrez Gene ID	283450
Protein Refseq	NP_001103132
UniProt ID	F8VWT9
Chromosome Location	12q24.13
Function	ligase activity; ubiquitin-protein transferase activity;