



Anti-CCDC57 polyclonal antibody (DPAB-DC1427)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC57 (coiled-coil domain containing 57) is a protein-coding gene.
Immunogen	Recombinant protein corresponding to amino acids of recombinant CCDC57. The sequence is HEEEISKDFPSSEIQRLREQNTSLRNAIAQMRKEMEALSHQIPPIQTAAESTDANQPDP EAGGDAATPDYVLALEAEIRTLKHKFKTLEKHLEDVLDPLKMSSPHAESQPSV
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC-P,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CCDC57 coiled-coil domain containing 57 [Homo sapiens (human)]
Official Symbol	CCDC57
Synonyms	CCDC57; coiled-coil domain containing 57; coiled-coil domain-containing protein 57;
Entrez Gene ID	284001
Protein Refseq	NP_932348
UniProt ID	Q2TAC2
Chromosome Location	17q25.3