



Anti-C1ORF204 polyclonal antibody (DPAB-DC1430)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C1orf204 (chromosome 1 open reading frame 204) is a protein-coding gene. An important paralog of this gene is VSIG2.
Immunogen	Recombinant protein corresponding to amino acids of recombinant C1orf204. The sequence is GLLVSSPHCEEPCAHSCAHPGLPPHLVHKLPLSYLQTQDTDAASRRINAPLAAGWSWLRLWLVT
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC-P,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C1orf204 chromosome 1 open reading frame 204 [Homo sapiens (human)]
Official Symbol	C1ORF204
Synonyms	C1ORF204; chromosome 1 open reading frame 204; uncharacterized protein LOC284677;
Entrez Gene ID	284677
Protein Refseq	NP_001127705
UniProt ID	Q5VU13
Chromosome Location	1q23.2