



# Anti-RABL3 (full length) polyclonal antibody (DPAB-DC1432)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | RABL3 (RAB, member of RAS oncogene family-like 3) is a protein-coding gene. GO annotations related to this gene include GTP binding.                                                                                                                                                                                                                        |
| <b>Immunogen</b>           | RABL3 (AAH20832, 1 a.a. ~ 236 a.a) full-length recombinant protein with GST tag. The sequence is<br>MASLDRVKVLVLGDSGVGKSSLVHLLCQNQVLGNPSWTVGCSVDVRVHDYKEGTPEEKTC<br>YIELWDVGGSVGSASSVKSTRAVFYNSVNGIIFVHDLTNKKSSQNLRRWSLEALNRDLVP<br>TGVLVTNGDYDQEQFADNQIPLLIGTKLDQIHETKRHEVLTTTAFLAEDFNPEEINLDC<br>TNPRYLAAGSSNAVKLSRFFDKVIEKRYFLREGNQIPGFPDRKRFGAGTLKSLHYD |
| <b>Source/Host</b>         | Mouse                                                                                                                                                                                                                                                                                                                                                       |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                                                                                                                                       |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                                                                                                                                                                                                |
| <b>Applications</b>        | WB (Recombinant protein), ELISA,                                                                                                                                                                                                                                                                                                                            |
| <b>Size</b>                | 50 µl                                                                                                                                                                                                                                                                                                                                                       |
| <b>Buffer</b>              | 50 % glycerol                                                                                                                                                                                                                                                                                                                                               |
| <b>Preservative</b>        | None                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage</b>             | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                    |

## GENE INFORMATION

**Gene Name** [RABL3 RAB, member of RAS oncogene family-like 3 \[ Homo sapiens \(human\) \]](#)

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Official Symbol</b>     | RABL3                                                                 |
| <b>Synonyms</b>            | RABL3; RAB, member of RAS oncogene family-like 3; rab-like protein 3; |
| <b>Entrez Gene ID</b>      | <a href="#">285282</a>                                                |
| <b>Protein Refseq</b>      | <a href="#">NP_776186</a>                                             |
| <b>UniProt ID</b>          | <a href="#">Q5HYI8</a>                                                |
| <b>Chromosome Location</b> | 3q13.33                                                               |
| <b>Function</b>            | GTP binding;                                                          |