



Anti-TMEM14A (full length) polyclonal antibody (DPAB-DC1461)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM14A (transmembrane protein 14A) is a protein-coding gene. Diseases associated with TMEM14A include juvenile myoclonic epilepsy, and idiopathic generalized epilepsy. An important paralog of this gene is TMEM14B.
Immunogen	TMEM14A (AAH15097, 1 a.a. ~ 99 a.a) full-length recombinant protein with GST tag. The sequence is MDLIGFGYAALVTFGSIFGYKRRGGVPSLIAGLFVGCLAGYGAYRVSNDKRDVKVSLFTA FFLATIMGVRFKRSKKIMPAGLVAGLSLMMILRLVLLLL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [TMEM14A transmembrane protein 14A \[Homo sapiens \(human\) \]](#)

Official Symbol	TMEM14A
Synonyms	TMEM14A; transmembrane protein 14A; PTD011; C6orf73;
Entrez Gene ID	28978
Protein Refseq	NP_054770
UniProt ID	Q9Y6G1
Chromosome Location	6p12.2