



Anti-COMMD5 (aa 125-224) polyclonal antibody (DPAB-DC1462)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	COMMD5 (COMM domain containing 5) is a protein-coding gene. Diseases associated with COMMD5 include hypertension. GO annotations related to this gene include protein binding.
Immunogen	COMMD5 (NP_054785, 125 a.a. ~ 224 a.a) partial recombinant protein with GST tag. The sequence is DLASVVFSGSRPLLDSPAQQGAWLPHVADFRWRVDVAISTALARSLQPSVLMQLKLSDD GSAYRFEVPTAKFQELRYSVALVLKEMADLEKRCERRLQD
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	COMMD5 COMM domain containing 5 [Homo sapiens (human)]
Official Symbol	COMMD5

Synonyms	COMMD5; COMM domain containing 5; HCARG; HT002; COMM domain-containing protein 5; hypertension-related calcium-regulated gene protein;
Entrez Gene ID	28991
Protein Refseq	NP_001074472
UniProt ID	Q9GZQ3
Chromosome Location	8q24.3
Function	protein binding;