



Anti-TMEM230 (full length) polyclonal antibody (DPAB-DC1467)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM230 (transmembrane protein 230) is a protein-coding gene. Diseases associated with TMEM230 include chronic lymphocytic leukemia, and leukemia.
Immunogen	C20orf30 (AAH09768, 1 a.a. ~ 120 a.a) full-length recombinant protein with GST tag. The sequence is MMPSRTNLATGIPSSKVKYSRLSSTDDGYIDLQFKKTPPKIPYKAIALATVLFLLIGAFLLI IIGSLLLSGYISKGGADRAVPVLIIGILVFLPGFYHLRIAYYASKGYRGYSYDDIPDFDD
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMEM230 transmembrane protein 230 [Homo sapiens (human)]
Official Symbol	TMEM230

Synonyms	TMEM230; transmembrane protein 230; HSPC274; C20orf30; dJ1116H23.2.1; UPF0414 transmembrane protein C20orf30;
Entrez Gene ID	29058
Protein Refseq	NP_001009923
UniProt ID	Q96A57
Chromosome Location	20p13