



Anti-TMOD3 (aa 280-352) polyclonal antibody (DPAB-DC1506)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMOD3 (tropomodulin 3 (ubiquitous)) is a protein-coding gene. Diseases associated with TMOD3 include estrogen excess, and transitional cell carcinoma. GO annotations related to this gene include actin binding and tropomyosin binding. An important paralog of this gene is TMOD1.
Immunogen	TMOD3 (NP_055362, 280 a.a. ~ 352 a.a) partial recombinant protein with GST tag. The sequence is RDNETLAELKIDNQRRQLGTAVELEMAKMLEENTNILKFGYQFTQQGPRTRAANAITKNN DLVRKRRVEGDHQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [TMOD3 tropomodulin 3 \(ubiquitous\) \[Homo sapiens \(human\) \]](#)

Official Symbol	TMOD3
Synonyms	TMOD3; tropomodulin 3 (ubiquitous); UTMOD; tropomodulin-3; ubiquitous tropomodulin;
Entrez Gene ID	29766
Protein Refseq	NP_055362
UniProt ID	Q9NYL9
Chromosome Location	15q21.1-q21.2
Function	actin binding; tropomyosin binding;