



Anti-NDRG2 polyclonal antibody (DPAB-DC1511)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NDRG2 (NDRG family member 2) is a protein-coding gene. Diseases associated with NDRG2 include motor peripheral neuropathy, and glioblastoma. GO annotations related to this gene include protein binding and molecular_function. An important paralog of this gene is NDRG4.
Specificity	Specific to recombinant protein GX2322. This antibody detects mNDRG2 protein.
Immunogen	Recombinant GST fusion protein corresponding to 137 amino acids of mouse NdrG2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Ndr2 N-myc downstream regulated gene 2 [Mus musculus (house mouse)]
Official Symbol	NDRG2
Synonyms	NDRG2; N-myc downstream regulated gene 2; Ndr2; SYLD; AI182517; AU040374; protein NDRG2; protein Ndr2; N-myc downstream regulated 2; N-myc downstream-regulated gene 2 protein;
Entrez Gene ID	29811
Protein Refseq	NP_001139431
UniProt ID	Q9QYG0
Chromosome Location	14 C2; 14