



Anti-ZFP292 polyclonal antibody (DPAB-DC1533)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX0006. This antibody detects mZFP292 protein.
Immunogen	Recombinant GST fusion protein corresponding to 136 mouse Zfp292. The sequence is PMGFEASFLKFLEESAVKQKKNSDRDHSNSGSKRGSHSSRRHVDKAAVAGSSHVCCKDSEIFVQFANPSKLQCSENVKIVLDKTLKDRSELVLKQLQEMKPTVSLKKLEVLSNNPDRTVLKEISIGKATGRGQY
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [Zfp292 zinc finger protein 292 \[Mus musculus \(house mouse\) \]](#)

Official Symbol	ZFP292
Synonyms	ZFP292; zinc finger protein 292; Zfp15; Zn-15; Zn-16; Zfp-15; Krox-10; AI449016; mKIAA0530; 5730450D02Rik; 5830493J20Rik; 9430062L07Rik; zinc finger protein 15;
Entrez Gene ID	30046
Protein Refseq	NP_038917
UniProt ID	Q9Z2U2
Chromosome Location	4 A5; 4 16.86 cM
Function	DNA binding; RNA polymerase II core promoter proximal region sequence-specific DNA binding transcription factor activity involved in positive regulation of transcription; metal ion binding; nucleic acid binding