



Anti-BPIFA4P (internal region) polyclonal antibody (DPAB-DC1583)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BPIFA4P (BPI fold containing family A, member 4, pseudogene) is a pseudogene, and is affiliated with the lncRNA class. Diseases associated with BPIFA4P include breast cancer.
Specificity	This antibody detects endogenous levels of total BASE protein.
Immunogen	A synthetic peptide corresponding to internal of human BASE.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P, IF, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	BPIFA4P BPI fold containing family A, member 4, pseudogene [Homo sapiens (human)]
------------------	---

Official Symbol	BPIFA4P
Synonyms	BPIFA4P; BPI fold containing family A, member 4, pseudogene; BASE;
Entrez Gene ID	317716
Chromosome Location	20q11.21