



Anti-6430573F11RIK polyclonal antibody (DPAB-DC1585)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX2099. This antibody detects mC8orf79 protein.
Immunogen	Recombinant GST fusion protein corresponding to 181 amino acids of mouse 6430573F11Rik.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	6430573F11Rik RIKEN cDNA 6430573F11 gene [Mus musculus (house mouse)]
Official Symbol	6430573F11RIK

Synonyms	6430573F11RIK; RIKEN cDNA 6430573F11 gene; Kiaa1456; putative methyltransferase KIAA1456;
Entrez Gene ID	319582
Protein Refseq	NP_795926
UniProt ID	Q80WQ4
Chromosome Location	8 A4; 8
Function	1-hydroxy-6-methoxypyrene methyltransferase activity; 1-hydroxypyrene methyltransferase activity; 1-phenanthrol methyltransferase activity; 2-phytyl-1,4-naphthoquinone methyltransferase activity