



Anti-ZBTB39 polyclonal antibody (DPAB-DC1588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZBTB39 (zinc finger and BTB domain containing 39) is a protein-coding gene. GO annotations related to this gene include DNA binding and zinc ion binding. An important paralog of this gene is ZBTB32.
Specificity	Specific to recombinant protein GX1511. This antibody detects mZBTB39 protein. It also recognizes human ZBTB39 protein.
Immunogen	Recombinant GST fusion protein corresponding to 161 amino acids of mouse Zbtb39.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Zbtb39 zinc finger and BTB domain containing 39 [Mus musculus (house mouse)]
Official Symbol	ZBTB39
Synonyms	ZBTB39; zinc finger and BTB domain containing 39; mKIAA0352; 7030401O21Rik; zinc finger and BTB domain-containing protein 39;
Entrez Gene ID	320080
Protein Refseq	NP_932152
UniProt ID	Q6PDK0
Chromosome Location	10 D3; 10
Function	metal ion binding; molecular_function;