



# Anti-Uricase (native) polyclonal antibody (DPAB-DC1589)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified Bacillus Uricase.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Bacillus
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

Gene Name	<a href="#">ABC3735 uricase [ Bacillus clausii KSM-K16 ]</a>
Official Symbol	ABC3735

<b>Synonyms</b>	ABC3735; uricase;
<b>Entrez Gene ID</b>	<a href="#">3202700</a>
<b>UniProt ID</b>	<a href="#">Q5WBJ3</a>
<b>Chromosome Location</b>	NC_006582.1 (3892824..3893849, complement)
<b>Pathway</b>	Microbial metabolism in diverse environments; Purine metabolism;