



Anti-ATG2A polyclonal antibody (DPAB-DC1614)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ATG2A (autophagy related 2A) is a protein-coding gene. Diseases associated with ATG2A include multiple endocrine neoplasia. GO annotations related to this gene include protein binding. An important paralog of this gene is ATG2B.
Specificity	Specific to recombinant protein GX0890. This antibody detects mATG2A protein.
Immunogen	Recombinant GST fusion protein corresponding to 213 amino acids of mouse Atg2a.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Atg2a autophagy related 2A [Mus musculus (house mouse)]
Official Symbol	ATG2A
Synonyms	ATG2A; autophagy related 2A; BC023754; mKIAA0404; A830054M12; 1810013C15Rik; autophagy-related protein 2 homolog A; ATG2 autophagy related 2 homolog A;
Entrez Gene ID	329015
Protein Refseq	NP_919329
UniProt ID	Q6P4T0
Chromosome Location	19 A; 19
Function	molecular_function;