



Anti-D630045J12RIK polyclonal antibody (DPAB-DC1620)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX0620. This antibody detects mKIAA1549 protein.
Immunogen	Recombinant GST fusion protein corresponding to 190 amino acids of mouse D630045J12Rik.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	D630045J12Rik RIKEN cDNA D630045J12 gene [Mus musculus (house mouse)]
Official Symbol	D630045J12RIK

Synonyms	D630045J12RIK; RIKEN cDNA D630045J12 gene; Kiaa1549; mKIAA1549; UPF0606 protein KIAA1549;
Entrez Gene ID	330286
Protein Refseq	NP_918950
UniProt ID	Q68FD9
Chromosome Location	6 B1; 6
Function	molecular_function;