



Anti-GLNA (native) polyclonal antibody (DPAB-DC1629)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Brevibacterium glnA.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Brevibacterium
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	glnA glutamine synthetase I [Corynebacterium glutamicum ATCC 13032]
Official Symbol	GLNA

Synonyms	GLNA; glutamine synthetase I;
Entrez Gene ID	3345165
UniProt ID	Q79VE3
Chromosome Location	NC_006958.1 (2320501..2321934)
Pathway	Alanine, aspartate and glutamate metabolism; Arginine and proline metabolism; Biosynthesis of amino acids; Glyoxylate and dicarboxylate metabolism