



Anti-ST7-OT4 polyclonal antibody (DPAB-DC1639)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ST7-OT4 (ST7 overlapping transcript 4) is a protein-coding gene.
Immunogen	Recombinant protein corresponding to amino acids of recombinant ST7OT4. The sequence is SQYCLLRFRQPCGAYSSCDEGKCVPGEVTESGAAFREQLQMAQASLPPQGSQRHHHNQAAR SQGFHRPSSVPVG
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ST7-OT4 ST7 overlapping transcript 4 [Homo sapiens (human)]
Official Symbol	ST7-OT4
Synonyms	ST7-OT4; ST7 overlapping transcript 4; ST7OT4; NCRNA00042;
Entrez Gene ID	338069
Chromosome Location	7q31.2