



Anti-C2ORF74 polyclonal antibody (DPAB-DC1651)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C2orf74 (chromosome 2 open reading frame 74) is a protein-coding gene. Diseases associated with C2orf74 include crohn's disease.
Immunogen	Recombinant protein corresponding to amino acids of recombinant C2orf74. The sequence is DANGGVDCAAAKVVTSNPEDHERILMQVMNLNVPMPRGILVQRQSKEVLATPLENRRDME AEEENQINEKQEPENAGETGQEEDDGLQKIHTSVTRTPSVVESQKRPLKGVTFSSREVIVV DLGNEYPTPRSYT
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C2orf74 chromosome 2 open reading frame 74 [Homo sapiens (human)]
Official Symbol	C2ORF74
Synonyms	C2ORF74; chromosome 2 open reading frame 74; uncharacterized protein C2orf74;
Entrez Gene ID	339804
Protein Refseq	NP_001137431
UniProt ID	A8MZ97
Chromosome Location	2p15