



Anti-TRIML1 (aa 369-468) polyclonal antibody (DPAB-DC1654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TRIML1 (tripartite motif family-like 1) is a protein-coding gene. GO annotations related to this gene include zinc ion binding and ligase activity. An important paralog of this gene is TRIML2.
Immunogen	FLJ36180 (NP_848651, 369 a.a. ~ 468 a.a) partial recombinant protein with GST tag. The sequence is PKPPGDLFSLIGLKIGDDYSLWSSPLKGQHVREPVCVKGVFLDYESGHIAFYNGTDESL IYSFPQASFQEALRPIFSPCLPNEGTTNDPLTICSLNSHV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TRIML1 tripartite motif family-like 1 [Homo sapiens (human)]
Official Symbol	TRIML1

Synonyms	TRIML1; tripartite motif family-like 1; RNF209; probable E3 ubiquitin-protein ligase TRIML1; RING finger protein 209; tripartite motif family-like protein 1;
Entrez Gene ID	339976
Protein Refseq	NP_848651
UniProt ID	Q8N9V2
Chromosome Location	4q35.2
Function	ligase activity; zinc ion binding;