



Anti-NACA2 (aa 116-215) polyclonal antibody (DPAB-DC1665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NACA2 (nascent polypeptide-associated complex alpha subunit 2) is a protein-coding gene. GO annotations related to this gene include TBP-class protein binding and transcription coactivator activity. An important paralog of this gene is NACAD.
Immunogen	NACAL (NP_954984, 116 a.a. ~ 215 a.a) partial recombinant protein with GST tag. The sequence is ASDAYIVFGGEAKIQDLSQQAQLAAAEKFRVQGEAVGNIQENTQTPTVQEESEEEVEDETG VEVKDVKLVMSQANVSRKAVRALKNNSNDIVNAIMELTV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NACA2 nascent polypeptide-associated complex alpha subunit 2 [Homo sapiens (human)]
------------------	---

Official Symbol	NACA2
Synonyms	NACA2; nascent polypeptide-associated complex alpha subunit 2; ANAC; NACAL; nascent polypeptide-associated complex subunit alpha-2; alpha-NAC protein; hom s 2.01 protein; IgE autoantigen alpha-nascent polypeptide-associated complex; nascent polypeptide-associated complex, alpha polypeptide, 2;
Entrez Gene ID	342538
Protein Refseq	NP_954984
UniProt ID	Q9H009
Chromosome Location	17q23.2